

Data Sheet

FluoTag-X2 anti-mEOS, Fluorophore-conjugated

Fluorophore	Catalog No. (20µl)	Catalog No. (200µl)
Abberior-Star 580	N3102-Ab580-S	N3102-Ab580-L
Abberior-Star RED	N3102-AbRED-S	N3102-AbRED-L
Abberior-Star 635P	N3102-Ab635P-S	N3102-Ab635P-L
AlexaFluor 647	N3102-AF647-S	N3102-AF647-L
Atto488	N3102-At488-S	N3102-At488-L
Atto542	N3102-At542-S	N3102-At542-L
Atto565	N3102-At565-S	N3102-At565-L
Atto647N	N3102-At647N-S	N3102-At647N-L
Sulfo-Cyanine 3	N3102-SC3-S	N3102-SC3-L
Sulfo-Cyanine 5	N3102-SC5-S	N3102-SC5-L

Identity	Camelid sdAb anti-mEOS; Clone 1E8
Content	2.5 µM fluorescently labeled sdAb in buffered saline, 50% glycerol, 0.09% sodium azide.
Storage	Up to 3 months: -20 °C Up to 12 months: -80 °C or below Protect from light!
Applications	IF/ICC Recommended dilution 1:500
Specificity	Recognizes mEOS and most common mEOS derivatives

For more information please visit our web page at www.nano-tag.com

To be used *in vitro* / for research only.
Non-toxic, non-hazardous, non-infectious.