

Data Sheet

FluoTag-X2 anti-GFAP Fluorophore-conjugated

Fluorophore	Catalog No. (200µl)
Alexa Fluor® 647	N3802-AF647-L
Abberior® STAR RED	N3802-AbRED-L
Abberior® STAR 580	N3802-Ab580-L
Abberior® STAR 635P	N3802-Ab635P-L
Atto488	N3802-At488-L
Atto542	N3802-At542-L
Atto565	N3802-At565-L
Atto647N	N3802-At647N-L
Sulfo-Cyanine 3	N3802-SC3-L
Sulfo-Cyanine 5	N3802-SC5-L

Identity	Camelid sdAb anti-GFAP, Clone 1D12
Content	2.5 µM labeled sdAb, 5 µM fluorophore in buffered saline, 50% glycerol, 0.09% sodium azide.
Storage	Up to 3 months: -20 °C Up to 12 months: -80 °C or below Protect from light!
Applications	IF/ICC/IHC Recommended dilution 1:500
Antigen	human GFAP
Specificity	Recognizes specifically GFAP from human, mouse and rat. Other species not tested.

For more information please visit our web page at www.nano-tag.com

To be used *in vitro* / for research only.
Non-toxic, non-hazardous, non-infectious.