
Data Sheet

SARS-CoV-2 Nucleocapsid protein (His-tagged)

Catalog no: P0101

Construct AA 1-419 from SARS-CoV-2 Nucleocapsid protein (UniProt ID: P0DTC9) fused to c-terminal polyhistidine tag

Content 1 mg lyophilized from 330 μ L phosphate buffer pH 6.0, 300 mM NaCl, 1 mM EDTA.
For reconstitution add 330 μ l H₂O, then aliquot and store at -20°C until use.

Storage Up to 6 months: -20 °C
Up to 12 months: -80 °C or below

Expressed in *E. coli*

Molecular Mass 46.4 kDa

Purity > 90 % based on SDS-PAGE analysis

By purchasing this product you agree to our general terms and conditions.

To be used *in vitro* / for research only.
Non-toxic, non-hazardous, non-infectious.