

FAQ - What is the concentration of antigen specific IgGs in polyclonal antiserum?

IgG polyclonal serum contains other proteins such as albumin and additional immunoglobulins which interfere with protein quantification. Consequently, only a minor part of the total IgG fraction is specific for the immunized antigen. A major part is directed against other unrelated antigens. Therefore, the exact concentration of antigen specific antibodies in polyclonal serum is unknown.

Service Unavailable

The server is temporarily unable to service your request due to maintenance downtime or capacity problems. Please try again later